

LIPOCONTROL

«EVERYTHING IS UNDER CONTROL»

The Latest Generation of Laser Lipolysis:

LipoControl is the new advanced device that allows optimization in the treatment of localized fat deposits. Thanks to its revolutionary technology of control and regulation, the laser system shows an accurate and controlled on-screen amount of energy delivery for maximized results


- Real-time visualization
- Mapping of the energy delivery
- Automatic adjustment of laser power
- Laser power off whenever the cannula stops moving



Results optimization: even distribution of the energy


Without control and regulation

(standard laser device)



With control and regulation

(LipoControl)




Onscreen mapping system shows real time distribution of energy

Advantages of Laser Lipolysis:

- A 5-year clinical insights
- 1 single session
- Local tumescent anesthesia
- No traumatic procedure
- Homogeneous results
- Skin tightening
- No social downtime

Areas:



Before / After:

BEFORE



AFTER



At 3 months. Photo courtesy of Dr Sadick



At 3 months. Photo courtesy of Dr Sadick



At 6 months. Photo courtesy of Dr Reynaud

Technical specifications:

Laser:	Diode laser
Wavelength:	980 nm
Power:	25 Watts
Emission mode:	Continuous, regulated
Laser Class:	Class 4
Dimensions:	180 x 70 x 64 cm
Weight:	55 kg

OSYRIS MEDICAL

Parc de la Haute Borne
60, Avenue Halley
59650 VILLENEUVE D'ASCQ - France

Tel.: +33 (0)3 20 67 90 00
Fax: +33 (0)3 20 04 46 24
Email: osyris.medical@osyris.com



www.osyrismedical.com

